

ADVISORY NOTICE

11_2024_UREA

Dear customers,

We inform you about the correction of a printing error in the IFUs of following products:

Cat.No.	Pack Name
XSYS0020	UREA 275
XSYS0075	UREA 564 XL - 1000
BLT00060	UREA 1000
BLT00061	UREA 250

The correction refers to changing the Limit of quantification and technical limit low in Assay parameters for XL analysers.

The current value 11.5 mg/dL (1.92 mol/L) is wrong.

The new correct value is 1.5 mg/dL (0.25 mol/L).

Please, change this value in your analyser.

In the case of ERBA XL automatic system contact your application specialist or service technician to update of assay parameters for UREA.

In case of any questions, please do not hesitate to contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.
Technical Product Management



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 21st May 2024

Ref. No. AN_11_2024