

ADVISORY NOTICE

16_2024_LAURA XL SW version 1.2.4.16

Dear customer,

We hereby inform you that the updated software (version 1.2.4.16) for LAURA XL automatic urine analyser was released on the 21st May. The updated SW version includes changes such as:

- new cleaning process customizable by user
- updated user manual integration
- perform system backup at each SW startup
- camera management improvements
- improvement on discrepancy on Colorimetry results for light color urines
- N/A flagged Results for Chemistry changed for more informative messages

and other changes based on market requests and bug corrections.

Erba Lachema highly recommends you update the SW on your instrument to the latest version (1.2.4.16).

For the latest SW version, please contact your local service engineer or distributor to install updated SW version to your instrument. The updated SW version you can find on ERBA Knowledge Base.

Applied for:

INSTRUMENT:

LAURA XL

cat.no.: INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management

V.Z.
Pavel

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno 3rd June 2024

Ref. No.AN 16_2024

