

ADVISORY NOTICE

22_2024_H360_SW_vA15.0

Dear customer,

We hereby inform you that the H360 device software has been updated to SW version A15.0 issued on 10th July 2024.

H360 SW_vA15.0 implements bugs fixes:

- Solved the problem where analyser could not input letter into the sample ID. SW_vA15.0 now supports input of numbers and letters to sample ID.
- Added a model name when printing quality control analysis and scattergram.
- SW_vA15.0 supports Ukrainian language.

SW version A15.0 is already included in devices produced from SN K10012422001 onwards.

For SW version A15.0, please contact your local service engineer or distributor to set up the new SW version from the web site [ERBA Knowledge Base](#).

Applied for:

INSTRUMENT:

H360

cat.no.INS00077

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-



Brno, 15th August 2024

Ref. No.AN_22_2024

