

## ADVISORY NOTICE

29\_2024\_XL-200 Plus\_SW\_v2024.02

Dear customer,

We hereby inform you that the XL-200 Plus device software has been updated to SW version 2024.02 issued on 9<sup>th</sup> May 2024.

XL-200 Plus SW\_v2024.02 includes a new feature:

- Preventive measure – Wash stirrer on COP (Carry over pair) wash.  
During carry over wash cycle stirrer will get wash in the respective cuvette/Trough.

SW version 2024.02 is already included in devices produced from SN 25P0124 onwards.

For SW version 2024.02, please contact your local service engineer or distributor.  
The latest version of the SW is available on the [ERBA Knowledge Base](#) website.


Applied for:

**INSTRUMENT:** XL-200 Plus cat.no : INS00001, INS000002

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Brno, 4<sup>th</sup> August 2024



Ref. No.AN\_29\_2024