

## ADVISORY NOTICE

28\_2024\_New sw version for LAURA XL – 1.2.5

Dear customer,

We hereby inform you that there was a new software version released for LAURA XL automatic urine analyser. This sw version is 1.2.5 and is available from 8<sup>th</sup> October 2024. This sw update includes enhancements to the system performance and bugfixes. Some of the updates include:

- Improvement of strips detection, allowing lower rate of rejected strips.
- Ergonomy improvement: QC printing title correction and multi selection for results deletion.

These enhancements are mainly focused for production of the instrument, therefore it is not compulsory to update to this software version.

Applied for:

**INSTRUMENT:**


LAURA XL

Cat. No:INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Brno, 14<sup>th</sup> October 2024

Ref. No.AN\_28\_2024

