

ADVISORY NOTICE

37_2024_Re-assignment of sediment for URINORM XL lot 2404127

Dear customer,

We hereby inform you that the **target values** for control material URINORM XL P (Cat. No. REG00060) lot **2404127** have been re-assigned.

The reassignment for evaluation of **RBCs in sediment** was carried out based on retesting of control material.

New target values for URINORM XL P are:

Lot	Expiration	URINORM XL P	
2404127	03/2026	2 – 30 / 11 – 158	N/hpf / N/ μ l

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

R&D specialist
Kateřina Seidlová**Erba Lachema s.r.o.**

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-

Brno, 13rd November

Ref. No.AN_37_2024