

ADVISORY NOTICE

38_2024_ New sw version for LAURA XL – 1.2.6.12

Dear customer,

We hereby inform you that there was a new software version released for LAURA XL automatic urine analyser. This sw version is 1.2.6.12 and is available from 20th January 2025. This sw update includes enhancements to the system performance and bugfixes, such as:

- focus detection improvement,
- improvement of wrong strip detection,
- improvement of colorimeter function,
- improvement in the cleaning of the transfer path,
- option for mandatory QC check before run was removed from the sw. QC measurements can be requested according to your internal Quality Control policy as before.

Applied for:

INSTRUMENT:

LAURA XL

Cat. No:INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 21st January 2025

Ref. No.AN_38_2024