

## ADVISORY NOTICE

02\_2025\_ECL 760\_SW\_v2.2e

Dear customer,

We hereby inform you that the device software of ECL 760 has been updated to SW version 2.2e issued on October 2024.

ECL 760 device SW update\_v2.2e implements the following changes:

The automatic calibration function of the optical path and its related functions are updated in the software, and the processing function of large data volume in the database is added.

SW version 2.2 e is already included in devices produced from SN 0439-71-291024 onwards.

For SW version 2.2e, please contact your local service engineer or distributor to set up the new SW version from the web site [ERBA Knowledge Base](#).

Applied for:

**INSTRUMENT:**

ECL 760

cat.no.INS00070

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

**Erba Lachema s.r.o.**

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-

Brno, February 3<sup>rd</sup>, 2025

Ref. No.AN\_02\_2025