

ADVISORY NOTICE

04_2025_ New SW version for LAURA

Dear customer,

We hereby inform you that the LAURA urine reader software has been updated to SW version 1.19 issued on 21st January 2025.

This SW update addresses a reported issue with the export of measured data to the LIS (Laboratory Information System).

The new SW version 1.19 is already included in devices produced from SN 2711 onwards.

The latest version of the SW is available on the [ERBA Knowledge Base](#) website.

Applied for:

INSTRUMENT:

LAURA

Cat. No: INS00063

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-



Brno, 6th February 2025

Ref. No.AN_04_2025